

Anti-Gstm5 Rabbit Polyclonal Antibody For Wb, Elisa - P46439

Item Number: CM0013152

Introduction

Anti-GSTM5 Rabbit Polyclonal Antibody validated for WB and ELISA, with cross reactivity to Mouse and Rat. Targets glutathione S-transferase Mu 5, a cytoplasmic enzyme involved in detoxification of electrophiles.

[Learn More](#)

Parameter	Value
Product Name	GSTM5 Rabbit pAb
Remarks/Alias	GTM5; GSTM5-5; GSTM5
Species	Human
GeneID (Human)	2949
GeneID	2949
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P46439
Protein Weight	26kDa
Shipping	Ice bag