

Anti-Gsto1 Polyclonal Antibody For Wb And Elisa - P78417

Item Number: CM0016756



Introduction

GSTO1 Rabbit pAb is a polyclonal antibody raised against recombinant human GSTO1 (aa 1-241). Suitable for WB and ELISA with cross-reactivity to human, mouse, rat. Ideal for arsenic metabolism and oxidative stress research.

[Learn More](#)

Parameter	Value
Product Name	GSTO1 Rabbit pAb
Remarks/Alias	P28; SPG-R; GSTO 1-1; GSTTLp28; HEL-S-21; GSTO1
Species	Human
GeneID	9446
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-241 of human GSTO1 (NP_004823.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P78417
Protein Weight	28kDa
Shipping	Ice bag