

Anti-Gsto2 Polyclonal Antibody For If/Icc And Elisa - Q9H4Y5

Item Number: CM0018431



Introduction

Rabbit polyclonal antibody against human GSTO2 (glutathione S-transferase omega-2), validated for IF/ICC and ELISA. Suitable for research on arsenic metabolism, ascorbate recycling, and oxidative stress.

[Learn More](#)

Parameter	Value
Product Name	GSTO2 Rabbit pAb
Remarks/Alias	GSTO 2-2; bA127L20.1; GSTO2
Species	Human
GeneID (Human)	119391
GeneID	119391
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 124-243 of human GSTO2 (NP_899062.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9H4Y5
Protein Weight	28kDa
Shipping	Ice bag