

Anti-Gsto2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9H4Y5

Item Number: CM0029525



Introduction

Rabbit polyclonal antibody against human GSTO2 (Glutathione S-transferase omega-2), validated for WB, IHC-P, IF/ICC, ELISA. Reactive with human, mouse, rat. Suitable for oxidative stress and arsenic metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	GSTO2 Rabbit pAb
Remarks/Alias	GSTO 2-2; bA127L20.1; GSTO2
Species	Human
GeneID (Human)	119391
GeneID	119391
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H4Y5
Protein Weight	28kDa
Shipping	Ice bag