


Anti-Gstt2 Rabbit Polyclonal Antibody For Wb And Elisa - P0CG29

Item Number: CM0014611



Introduction

Rabbit polyclonal antibody targeting human GSTT2 (Glutathione S-transferase theta-2). Validated for WB and ELISA, cross-reacts with human, mouse, and rat. Ideal for detoxification and cancer metabolism research.

[Learn More](#)

Parameter	Value
Product Name	GSTT2 Rabbit pAb
Remarks/Alias	GSTT2
Species	Human
GeneID (Human)	2953
GeneID	2953
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-244 of human GSTT2 (NP_000845.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P0CG29
Protein Weight	28kDa
Shipping	Ice bag