


Anti-Gstt2B Polyclonal Antibody For Wb, If/Icc, Elisa - P0Cg30

Item Number: CM0016906



Introduction

High-quality rabbit polyclonal antibody against GSTT2B (GSTT2P). Validated for WB, IF/ICC, ELISA with cross-reactivity to human, mouse, rat. Ideal for studying glutathione conjugation and detoxification pathways.

[Learn More](#)

Parameter	Value
Product Name	GSTT2B Rabbit pAb
Remarks/Alias	GSTT2P; GSTT2B
Species	Human
Gene ID (Human)	653689
Gene ID	653689
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 115-244 of human GSTT2B (NP_001074312.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P0CG30
Protein Weight	28kDa
Shipping	Ice bag