


Anti-Gtf2A1 Rabbit Polyclonal Antibody For Wb And Elisa - P52655

Item Number: CM0025831



Introduction

Rabbit polyclonal antibody targeting human GTF2A1 (TFIIA subunit 1) for Western blot and ELISA. Detects a 42kDa nuclear protein involved in RNA polymerase II transcription initiation complex.

[Learn More](#)

Parameter	Value
Product Name	GTF2A1 Rabbit pAb
Remarks/Alias	TF2A1; TFIIA; TFIIAL; TFIIA-42; GTF2A1
Species	Human
GeneID (Human)	2957
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human GTF2A1 (NP_056943.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P52655
Protein Weight	42kDa
Shipping	Ice bag