


Anti-Gtf2A2 Rabbit Polyclonal Antibody For Wb And Elisa - P52657

Item Number: CM0028008



Introduction

Rabbit polyclonal antibody targeting human GTF2A2 (TFIIA-gamma, TFIIA-12, TFIIAS). Confirmed for Western blot and ELISA applications. Ideal for RNA polymerase II transcription research.

[Learn More](#)

Parameter	Value
Product Name	GTF2A2 Rabbit pAb
Remarks/Alias	TF2A2; TFIIA; T18745; TFIIAS; HsT18745; TFIIA-12; TFIIA-gamma; GTF2A2
Species	Human
GeneID (Human)	2958
GeneID	2958
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-109 of human GTF2A2 (NP_004483.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P52657
Protein Weight	12kDa
Shipping	Ice bag