


Anti-Gtf2E1 Polyclonal Antibody For Wb, If, Elisa - P29083

Item Number: CM0016649



Introduction

Rabbit polyclonal antibody against human GTF2E1, validated for Western blot, immunofluorescence/ICC, and ELISA. Cross-reacts with mouse and rat. Recombinant immunogen, approximately 49 kDa. Suitable for transcription research.

[Learn More](#)

Parameter	Value
Product Name	GTF2E1 Rabbit pAb
Remarks/Alias	FE; TF2E1; TFIIE-A; GTF2E1
Species	Human
GeneID (Human)	2960
GeneID	2960
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 230-439 of human GTF2E1 (NP_005504.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29083
Protein Weight	49kDa
Shipping	Ice bag