


Anti-Gtf2E2 Polyclonal Antibody For Wb, Elisa - P29084

Item Number: CM0027163



Introduction

This rabbit polyclonal antibody targets human transcription initiation factor IIE subunit beta (GTF2E2). Suitable for WB and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	GTF2E2 Rabbit pAb
Remarks/Alias	FE; TTD6; TF2E2; TFIIE-B; GTF2E2
Species	Human
Gene ID (Human)	2961
Gene ID	2961
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-291 of human GTF2E2 (NP_002086.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29084
Protein Weight	33kDa
Shipping	Ice bag