


Anti-Gtf2F2 Polyclonal Antibody For Wb, If/Icc, Elisa - P13984

Item Number: CM0025180



Introduction

Rabbit polyclonal antibody targeting human GTF2F2 (TFIIF-beta), validated for WB, IF/ICC, and ELISA. Suitable for studying general transcription initiation factor complexes.

[Learn More](#)

Parameter	Value
Product Name	GTF2F2 Rabbit pAb
Remarks/Alias	BTF4, RAP30, TF2F2, TFIIF, GTF2F2
Species	Human
Gene ID (Human)	2963
Gene ID	2963
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-249 of human GTF2F2 (NP_004119.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P13984
Protein Weight	28kDa
Shipping	Ice bag