


Anti-Gtf2f2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P13984

Item Number: CM0002371



Introduction

Rabbit monoclonal antibody against human GTF2F2 (TFIIF-beta, RAP30) validated for Western blot, immunohistochemistry, and ELISA. Reacts with human, mouse, and rat. For transcription initiation research.

[Learn More](#)

Parameter	Value
Product Name	GTF2F2 Rabbit mAb
Remarks/Alias	BTf4; RAP30; TF2F2; TFIIF; GTF2F2
Species	Human
GeneID (human)	2963
GeneID	2963
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GTF2F2 (P13984).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13984
Protein Weight	28kDa
Shipping	Ice bag