


Anti-Gtf2H1 Polyclonal Antibody For Wb, If/Icc, Elisa - P32780

Item Number: CM0027620



Introduction

Rabbit polyclonal antibody against human GTF2H1 (p62, BTF2) for Western blot, immunofluorescence/immunocytochemistry, and ELISA. Validated in human, mouse, and rat samples. Suitable for transcription and DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	GTF2H1 Rabbit pAb
Remarks/Alias	P62; BTF2; TFB1; TFIIH; GTF2H1
Species	Human
Gene ID (Human)	2965
Gene ID	2965
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 1-425 of human GTF2H1 (NP_001135779.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P32780
Protein Weight	62kDa
Shipping	Ice bag