


Anti-Gtf2H3 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q13889

Item Number: CM0028737



Introduction

Rabbit polyclonal antibody targeting human GTF2H3, the TFIIH core complex subunit (34kDa). Validated for WB, IP, and ELISA; reacts with human and mouse. Ideal for studying transcription initiation and nucleotide excision repair mechanisms.

[Learn More](#)

Parameter	Description
Product Name	GTF2H3 Rabbit pAb
Remarks/Alias	P34; BTF2; TFB4; TFIIH; GTF2H3
Species	Human
Gene ID	2967
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-308 of human GTF2H3 (NP_001507.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13889
Protein Weight	34kDa
Shipping	Ice bag