


Anti-Gtf2H4 Rabbit Polyclonal Antibody For Wb, Elisa - Q92759

Item Number: CM0027624



Introduction

Rabbit polyclonal antibody targeting human GTF2H4 (TFIIH p52 subunit), suitable for Western blot and ELISA, cross-reactive with human and mouse. A valuable tool for transcription and DNA repair research.

[Learn More](#)

Parameter	Specification
Product Name	GTF2H4 Rabbit pAb
Remarks/Alias	P52; TFB2; TFIIH; GTF2H4
Species	Human
GeneID (Human)	2968
GeneID	2968
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human GTF2H4 (NP_001508.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92759
Protein Weight	52kDa
Shipping	Ice bag