

Anti-Gtf2Ird1 Polyclonal Antibody For Wb, Elisa - Q9Uhl9

Item Number: CM0027210



Introduction

Rabbit polyclonal antibody against GTF2IRD1, suitable for Western blot and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for transcription regulation and muscle differentiation research.

[Learn More](#)

Parameter	Value
Product Name	GTF2IRD1 Rabbit pAb
Remarks/Alias	BEN; WBS; GTF3; RBAP2; CREAM1; MUSTRD1; WBSCR11; WBSCR12; hMusTRD1alpha1; GTF2IRD1
Species	Human
Gene ID (Human)	9569
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 660-959 of human GTF2IRD1 (NP_057412.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UHL9
Protein Weight	106kDa
Shipping	Ice bag