

Anti-Gtf2Ird1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Uhl9

Item Number: CM0021707



Introduction

High-quality rabbit polyclonal antibody against human GTF2IRD1, validated for Western blot and ELISA. Cross-reacts with rat. Ideal for transcription regulation, muscle differentiation, and Williams-Beuren syndrome research.

[Learn More](#)

Parameter	Value
Product Name	GTF2IRD1 Rabbit pAb
Remarks/Alias	GTF2IRD1; BEN; CREAM1; GTF3; MUSTRD1; RBAP2; WBS; WBSER11; WBSER12; hMusTRD1alpha1; GTF2I repeat domain containing 1
Species	Human
GeneID	9569
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 660-890 of human GTF2IRD1 (NP_057412.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q9UHL9
Protein Weight	104kDa/106kDa/107kDa
Shipping	Ice bag