

Anti-Gtf3C3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc - Q9Y5Q9

Item Number: CM0020459



Introduction

Rabbit polyclonal antibody targeting human GTF3C3 (TFIIIC102, TFIIICγ). Validated for WB, IHC-P, IF/ICC, ELISA. Detects mouse and rat GTF3C3. Suitable for RNA polymerase III transcription studies.

[Learn More](#)

Parameter	Value
Product Name	GTF3C3 Rabbit pAb
Remarks/Alias	TFIIIC102; TFIIICgamma; TFIIIC2-102; GTF3C3
Species	Human
GeneID (Human)	9330
GeneID	9330
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human GTF3C3 (NP_036218.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5Q9
Protein Weight	101kDa
Shipping	Ice bag