

Anti-Gtf3C4 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Ukn8

Item Number: CM0025524



Introduction

Rabbit polyclonal antibody against human GTF3C4 (TFIIIC90) for Western blot, immunofluorescence, and ELISA. Validated in human and rat. Essential for RNA polymerase III transcription and histone acetylation.

[Learn More](#)

Parameter	Specification
Product Name	GTF3C4 Rabbit pAb
Remarks/Alias	KAT12; TFIIIC90; TFIIIC90; TFIIIC290; TF3C-delta; TFIIICDELTA; GTF3C4
Species	Human
GeneID (human)	9329
GeneID	9329
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 583-822 of human GTF3C4 (NP_036336.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9UKN8
Protein Weight	92kDa
Shipping	Ice bag