


Anti-Gtpbp1 Rabbit Polyclonal Antibody For Wb, Elisa - 000178

Item Number: CM0017807



Introduction

Rabbit polyclonal antibody against human GTPBP1, validated for WB and ELISA. Cross-reacts with mouse. Targets a key translation elongation factor and RNA quality control protein.

[Learn More](#)

Parameter	Value
Product Name	GTPBP1 Rabbit pAb
Remarks/Alias	GP1, GP-1, HSPC018, GTPBP1
Species	Human
GeneID (Human)	9567
GeneID	9567
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 10-280 of human GTPBP1 (NP_004277.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	000178
Protein Weight	72kDa
Shipping	Ice bag