


Anti-Gtpbp10 Rabbit Polyclonal Antibody For Wb And Elisa - A4D1E9

Item Number: CM0019336



Introduction

GTPBP10 Rabbit pAb validated for WB and ELISA. Targets human GTPBP10 (ObgH2), a nucleolar GTPase essential for ribosome maturation. Ideal for studies on nucleolar architecture and ribosome biogenesis.

[Learn More](#)

Parameter	Value
Product Name	GTPBP10 Rabbit pAb
Remarks/Alias	ObgH2, UG0751c10, GTPBP10
Species	Human
GeneID	85865
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-387 of human GTPBP10 (NP_149098.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	A4D1E9
Protein Weight	43kDa
Shipping	Ice bag