


Anti-Gtpbp4 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Bze4

Item Number: CM0008328



Introduction

GTPBP4 Rabbit monoclonal antibody validated for WB, IHC-P, ELISA. Recognizes human GTPBP4, a nucleolar GTPase involved in 60S ribosomal biogenesis and TP53 regulation. Suitable for cancer and ribosomal research.

[Learn More](#)

Parameter	Value
Product Name	GTPBP4 Rabbit mAb
Remarks/Alias	NGB; CRFG; NOG1; GTPBP4
Species	Human
GeneID	23560
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2-92 of human GTPBP4 (Q9BZE4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS Accession	Q9BZE4
Protein Weight	74kDa
Shipping	Ice bag