


Anti-Gtpbp6 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O43824

Item Number: CM0021857



Introduction

GTPBP6 Rabbit pAb is a polyclonal antibody validated for WB, IF/ICC, and ELISA, cross-reacting with human, mouse, and rat. Ideal for studying putative GTP-binding protein 6 expression.

[Learn More](#)

Parameter	Value
Product Name	GTPBP6 Rabbit pAb
Remarks/Alias	PGPL; GTPBP6
Species	Human
GeneID (Human)	8225
GeneID	8225
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 140-350 of human GTPBP6 (NP_036359.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43824
Protein Weight	57kDa
Shipping	Ice bag