


Anti-Gtse1 Polyclonal Antibody For Wb And Elisa - Q9Nyz3

Item Number: CM0017444



Introduction

High-quality rabbit polyclonal antibody against human GTSE1, validated for WB and ELISA, with cross-reactivity to mouse and rat. Targets a cell cycle regulator involved in G2/M arrest.

[Learn More](#)

Parameter	Value
Product Name	GTSE1 Rabbit pAb
Remarks/Alias	B99; GTSE1
Species	Human
GeneID (Human)	51512
GeneID	51512
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-290 of human GTSE1 (NP_057510.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NYZ3
Protein Weight	77kDa
Shipping	Ice bag