

Anti-Gtsf1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Ww33

Item Number: CM0019161

Introduction

GTSF1 rabbit polyclonal antibody for WB, IF/ICC, ELISA. Validated in human, mouse, rat. Targets gametocyte-specific factor 1, essential for spermatogenesis and retrotransposon silencing. Predicted 19kDa molecular weight.

[Learn More](#)

Parameter	Value
Product Name	GTSF1 Rabbit pAb
Remarks/Alias	Cue110; FAM112B; GTSF1
Species	Human
GeneID (Human)	121355
GeneID	121355
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-167 of human GTSF1 (NP_653195.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WW33
Protein Weight	19kDa
Shipping	Ice bag