

Anti-Guanine Deaminase (Gda) Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y2T3

Item Number: CM0030360



Introduction

Rabbit polyclonal antibody against human GDA (guanine deaminase), suitable for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studying guanine deamination pathways in purine metabolism research.

[Learn More](#)

Parameter	Value
Product Name	GDA Rabbit pAb
Remarks/Alias	GAH; CYPIN; GUANASE; NEDASIN; GDA
Species	Human
GeneID	9615
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 165-454 of human GDA (NP_004284.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y2T3
Protein Weight	51kDa
Shipping	Ice bag