

Anti-Guanylyl Cyclase Beta 1 (Gucy1B3) Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q02153

Item Number: CM0006036



Introduction

Rabbit monoclonal antibody against Guanylyl Cyclase beta 1 (GUCY1B3), validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for nitric oxide-cGMP signaling research.

[Learn More](#)

Parameter	Value
Product Name	Guanylyl Cyclase beta 1 (GUCY1B3) Rabbit mAb
Remarks/Alias	GUCB3; GC-SB3; GUC1B3; GUCSB3; GUCY1B3; GC-S-beta-1; Guanylyl Cyclase beta 1 (GUCY1B3)
Species	Human
GeneID (Human)	2983
GeneID	2983
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02153
Protein Weight	71kDa
Shipping	Ice bag