

Anti-Guanylyl Cyclase Beta 1 (Gucy1B3) Rabbit Polyclonal Antibody For Wb And Elisa - Q02153

Item Number: CM0013653



Introduction

Rabbit polyclonal antibody against human Guanylyl Cyclase beta 1 (GUCY1B3), suitable for Western blot and ELISA. Cross-reacts with mouse and rat. Detects the nitric oxide receptor cGMP synthase subunit. Available from CamelBio, your IVD raw material partner.

[Learn More](#)

Parameter	Value
Product Name	Guanylyl Cyclase beta 1 (GUCY1B3) Rabbit pAb
Remarks/Alias	GUCB3; GC-SB3; GUC1B3; GUC5B3; GUCY1B3; GC-S-beta-1; Guanylyl Cyclase beta 1 (GUCY1B3)
Species	Human
Gene ID (Human)	2983
Gene ID	2983
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human Guanylyl Cyclase beta 1 (Guanylyl Cyclase beta 1 (GUCY1B3)) (NP_000848.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02153
Protein Weight	71kDa
Shipping	Ice bag