


Anti-Guca2A Polyclonal Antibody For Wb And Elisa - Q02747

Item Number: CM0027802



Introduction

Rabbit polyclonal antibody against human guanylin (GUCA2A), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Suitable for gastrointestinal research and enterotoxin signaling studies.

[Learn More](#)

Parameter	Value
Product Name	GUCA2A Rabbit pAb
Remarks/Alias	GUCA2; STARA; GCAP-I; GUCA2A
Species	Human
GeneID (Human)	2980
GeneID	2980
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-115 of human GUCA2A (NP_291031.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02747
Protein Weight	12kDa
Shipping	Ice bag