


Anti-Guca2B Polyclonal Antibody For Wb, Elisa - Q16661

Item Number: CM0027618



Introduction

Rabbit polyclonal antibody against human GUCA2B (uroguanylin, GCAP-II). Validated for WB and ELISA, cross-reacts with mouse. Useful for studying intestinal guanylate cyclase signaling and fluid/electrolyte transport.

[Learn More](#)

Parameter	Value
Product Name	GUCA2B Rabbit pAb
Remarks/Alias	UGN; GCAP-II; GUCA2B
Species	Human
GeneID (Human)	2981
GeneID	2981
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GUCA2B (NP_009033.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q16661
Protein Weight	12kDa
Shipping	Ice bag