

Anti-Gucy1A3 Polyclonal Antibody For Wb, If/Icc, Elisa - Q02108

Item Number: CM0014466



Introduction

Rabbit polyclonal antibody targeting human GUCY1A3, validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for studies of soluble guanylate cyclase signaling.

[Learn More](#)

Parameter	Value
Product Name	GUCY1A3 Rabbit pAb
Remarks/Alias	GUCA3; MYMY6; GC-SA3; GUC1A3; GUCSA3; GUCY1A3; GCS-alpha-3; GC-S-alpha-1
Species	Human
Gene ID	2982
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human GUCY1A3 (NP_000847.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02108
Protein Weight	77kDa
Shipping	Ice bag