

Anti-Gucy2D Polyclonal Antibody For Wb And Elisa - Q02846

Item Number: CM0028332

Introduction

High-quality rabbit polyclonal antibody against retinal guanylyl cyclase 1 (GUCY2D), suitable for Western blot and ELISA applications.

Demonstrates cross-reactivity with human, mouse, and rat samples. Supports vision and phototransduction studies.

[Learn More](#)

Parameter	Specification
Product Name	GUCY2D Rabbit pAb
Remarks/Alias	LCA; CACD; CG-E; CYGD; LCA1; RCD2; CACD1; CORD5; CORD6; GUC2D; ROSGC; retGC; CSNB1i; GUC1A4; RETGC-1; ROS-GC1; GUCY2D
Species	Human
Gene ID (Human)	3000
Gene ID	3000
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 52-350 of human GUCY2D (NP_000171.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02846
Protein Weight	120kDa
Shipping	Ice bag