

Anti-Gucy2F Rabbit Polyclonal Antibody For Wb And Elisa - P51841

Item Number: CM0018418

Introduction

Rabbit polyclonal antibody against GUCY2F (Retinal guanylyl cyclase 2), suitable for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studying cGMP synthesis in photoreceptor cells.

[Learn More](#)

Parameter	Value
Product Name	Anti-GUCY2F Rabbit Polyclonal Antibody
Remarks/Alias	CYGF; GC-F; GUC2F; GUC2DL; RETGC-2; ROS-GC2; GUCY2F
Species	Human
Gene ID (Human)	2986
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 250-460 of human GUCY2F (NP_001513.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P51841
Protein Weight	125kDa
Shipping	Ice bag