

Anti-Guk1 Rabbit Polyclonal Antibody For Wb, Elisa - Q16774

Item Number: CM0031009



Introduction

Guanylate kinase (GUK1) rabbit polyclonal antibody validated for Western blot and ELISA. Targets human GUK1; recombinant immunogen. Suitable for studying nucleotide metabolism, cGMP recycling, and prodrug activation.

[Learn More](#)

Parameter	Value
Product Name	GUK1 Rabbit pAb
Remarks/Alias	GMK; GUK1
Species	Human
GeneID	2987
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-197 of human GUK1 (NP_000849.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q16774
Protein Weight	22kDa
Shipping	Ice bag