

Anti-Gulp1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Ubp9

Item Number: CM0029370



Introduction

Rabbit polyclonal antibody targeting GULP1 (CED6, GULP) for WB, IHC-P, ELISA applications. Cross-reacts with human, mouse, and rat. Specific for phagocytosis, apoptosis, and lipid transport studies.

[Learn More](#)

Parameter	Value
Product Name	GULP1 Rabbit pAb
Remarks/Alias	CED6; GULP; CED-6; GULP1
Species	Human
GeneID (human)	51454
GeneID	51454
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human GULP1 (NP_057399.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBP9
Protein Weight	34kDa
Shipping	Ice bag