

Anti-Gw5 Polyclonal Antibody For Wb, Elisa - A0A1D8Gzc0

Item Number: CM0022665

Introduction

Rabbit polyclonal anti-GW5 antibody validated for Western blot and ELISA. Detects Oryza sativa GW5 calmodulin-binding protein (~49 kDa). Ideal for rice research and agricultural diagnostic applications.

[Learn More](#)

Parameter	Value
Product Name	GW5 Rabbit pAb
Remarks/Alias	GW5
Species	Oryza sativa
GeneID	4338011
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-300 of Oryza sativa GW5 (A0A1D8GZC0)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Oryza sativa
SWISS	A0A1D8GZC0
Protein Weight	49kDa
Shipping	Ice bag