


Anti-Gypb Polyclonal Antibody For If/Icc, Elisa - P06028

Item Number: CM0016083



Introduction

This rabbit polyclonal antibody targets human Glycophorin-B (GYPB/CD235b), a key erythrocyte membrane sialoglycoprotein. Validated for immunofluorescence (IF/ICC) and ELISA, supporting research on red blood cell membrane stability and blood group antigens.

[Learn More](#)

Parameter	Value
Product Name	GYPB Rabbit pAb
Remarks/Alias	SS; GPB; GYP; MNS; GYPA; PAS-3; CD235b; GYPB
Species	Human
GeneID (Human)	2994
GeneID	2994
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-91 of human GYPB (NP_002091.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P06028
Protein Weight	10kDa
Shipping	Ice bag