


Anti-Gypb Rabbit Polyclonal Antibody For Wb, Elisa - P06028

Item Number: CM0028424



Introduction

Rabbit polyclonal antibody against human Glycophorin-B (GYPB, CD235b), suitable for WB and ELISA. Cross-reacts with human and mouse. Ideal for erythrocyte membrane studies.

[Learn More](#)

Parameter	Value
Product Name	GYPB Rabbit pAb
Remarks/Alias	SS, GPB, GYP, MNS, GYPA, PAS-3, CD235b, GYPB
Species	Human
Gene ID (Human)	2994
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-91 of human GYPB (NP_002091.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P06028
Protein Weight	10kDa
Shipping	Ice bag