

Anti-Gzma Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, Elisa - P12544

Item Number: CM0012447



Introduction

High-quality rabbit monoclonal antibody against human Granzyme A (GZMA). Validated for WB, IP, IF/ICC, and ELISA. Suitable for studies on cytotoxic T-cell and NK-cell mediated pyroptosis and target cell death.

[Learn More](#)

Parameter	Value
Product Name	GZMA Rabbit mAb
Remarks/Alias	HFSP; CTLA3
Species	Human
GeneID (Human)	3001
GeneID	3001
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-262 of human GZMA (NP_006135.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P12544
Protein Weight	29kDa
Shipping	Ice bag