


Anti-Gzma Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P12544

Item Number: CM0028684



Introduction

Rabbit polyclonal antibody against human GZMA, validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Targets granzyme A, a serine protease involved in pyroptosis and CTL/NK cell cytotoxicity.

[Learn More](#)

Parameter	Value
Product Name	GZMA Rabbit pAb
Remarks/Alias	HFSP; CTLA3; GZMA
Species	Human
GeneID (Human)	3001
GeneID	3001
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-262 of human GZMA (NP_006135.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12544
Protein Weight	29kDa
Shipping	Ice bag