

Anti-H2-Ab1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P06345

Item Number: CM0021053

Introduction

Rabbit polyclonal antibody targeting mouse H2-Ab1 (MHC class II beta chain). Validated for WB, IHC-P, IF/ICC, and ELISA applications. Cross-reacts with human, mouse, and rat. Ideal for immunology and histocompatibility research.

[Learn More](#)

Parameter	Value
Product Name	H2-Ab1 Rabbit pAb
Remarks/Alias	IAb; Ia2; Ia-2; Abeta; H-2Ab; H2-Ab; Rmcs1; I-Abeta; H2-Ab1
Species	Mouse
GeneID (human)	14961
GeneID	14961
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-219 of mouse H2-Ab1 (NP_996988.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06345
Protein Weight	30kDa
Shipping	Ice bag