


Anti-H3F3B Rabbit Polyclonal Antibody For Wb, Elisa - P84243

Item Number: CM0014965



Introduction

H3F3B Rabbit Polyclonal Antibody is designed for WB and ELISA applications, specifically recognizing human Histone H3.3. This antibody targets the variant histone H3.3, ideal for chromatin and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	H3F3B Rabbit pAb
Remarks/Alias	H3-3A; H3.3B; H3F3B; BRYLIB2
Species	Human
Gene ID	3021
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-136 of human H3F3B (NP_002098.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P84243
Protein Weight	15kDa
Shipping	Ice bag