

Anti-H3Q5Ser (Histone H3.1, P68431) Rabbit Polyclonal Antibody For Wb, If, Elisa

Item Number: CM0000083



Introduction

Rabbit polyclonal antibody targeting serotonylated histone H3.1 at Q5 (P68431). Validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for epigenetic research, chromatin studies, and histone modification analysis.

[Learn More](#)

Parameter	Value
Product Name	H3Q5ser Rabbit pAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; H3Q5ser
Species	Human
GeneID (human)	8350
GeneID	8350
Immunogen	A synthetic serotonylated peptide around Q5 of human histone H3 (NP_003520.1)
Source	Rabbit
Category	Other Modified Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P68431
Protein Weight	15 kDa
Shipping	Ice bag