

Anti-Hace1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Iyu2

Item Number: CM0027412



Introduction

Rabbit polyclonal antibody targeting HACE1, validated in WB, IHC-P, and ELISA. Recognizes human, mouse, and rat HACE1 (102 kDa). E3 ubiquitin ligase involved in Golgi membrane dynamics, RAC1 degradation, and host defense. Ideal for cell biology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	HACE1 Rabbit pAb
Remarks/Alias	SPPRS; HACE1
Species	Human
GeneID (Human)	57531
GeneID	57531
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HACE1 (NP_065822.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IYU2
Protein Weight	102kDa
Shipping	Ice bag