

Anti-Hadh Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q16836

Item Number: CM0013296



Introduction

High-quality rabbit polyclonal antibody against human HADH, validated for western blot, IHC-paraffin, and ELISA. Cross-reacts with mouse and rat. Ideal for studying mitochondrial fatty acid beta-oxidation and metabolic disorders.

[Learn More](#)

Parameter	Value
Product Name	HADH Rabbit pAb
Remarks/Alias	HAD; HCDH; HHF4; HADH1; HADHSC; MSCHAD; HADH
Species	Human
GeneID (Human)	3033
GeneID	3033
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-314 of human HADH (NP_005318.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16836
Protein Weight	34kDa
Shipping	Ice bag