

Anti-Hadha Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P40939

Item Number: CM0029301



Introduction

This rabbit polyclonal antibody specifically recognizes HADHA (TP-alpha), validated for WB, IF/ICC, IP, ELISA, and cross-reacts with human, mouse, and rat. Ideal for studying mitochondrial beta-oxidation and fatty acid metabolism. ~83 kDa.

[Learn More](#)

Parameter	Specification
Product Name	[KO Validated] HADHA Rabbit pAb
Remarks/Alias	GBP; ECHA; HADH; LCEH; MTPA; LCHAD; TP-ALPHA; HA
Species	Human
GeneID	3030
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 484-763 of human HADHA (NP_000173.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40939
Protein Weight	83kDa
Shipping	Ice bag