

Anti-Hadha Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P40939

Item Number: CM0023346



Introduction

Rabbit polyclonal HADHA antibody validated for WB, IHC-P, IF/ICC, and ELISA. Detects human and mouse HADHA, an 83 kDa mitochondrial trifunctional enzyme subunit involved in fatty acid beta-oxidation.

[Learn More](#)

Parameter	Value
Product Name	HADHA Rabbit pAb
Remarks/Alias	GBP; ECHA; HADH; LCEH; MTPA; LCHAD; TP-ALPHA; HADHA
Species	Human
GeneID (Human)	3030
GeneID	3030
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 545-763 of human HADHA (NP_000173.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P40939
Protein Weight	83kDa
Shipping	Ice bag