

Anti-Hadha Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P40939

Item Number: CM0007375



Introduction

Rabbit monoclonal antibody targeting human HADHA, validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Detects mitochondrial trifunctional enzyme subunit alpha for research use.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] HADHA Rabbit mAb
Remarks/Alias	GBP; ECHA; HADH; LCEH; MTPA; LCHAD; TP-ALPHA; HA
Species	Human
Gene ID (Human)	3030
Gene ID	3030
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 545-763 of human HADHA (NP_000173.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40939
Protein Weight	83kDa
Shipping	Ice bag