

Anti-Hadhb Monoclonal Antibody (Rabbit) For Wb, Ihc-P, If/Icc, Elisa - P55084

Item Number: CM0012724



Introduction

High-quality rabbit monoclonal antibody against HADHB (Trifunctional enzyme subunit beta, mitochondrial) for WB, IHC-P, IF/ICC, and ELISA. Validated in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] HADHB Rabbit mAb
Remarks/Alias	ECHB; MTPB; MTPD; MTPD2; MSTP029; TP-BETA; HADHB
Species	Human
Gene ID (Human)	3032
Gene ID	3032
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 331-474 of human HADHB (NP_000174.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55084
Protein Weight	51kDa
Shipping	Ice bag