


Anti-Hagh Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q16775

Item Number: CM0025915



Introduction

High-quality HAGH rabbit polyclonal antibody validated for western blot, immunofluorescence, and ELISA. Detects human, mouse, and rat hydroxyacylglutathione hydrolase. Ideal for studying glyoxalase II in mitochondrial and cytosolic detoxification pathways.

[Learn More](#)

Parameter	Value
Product Name	HAGH Rabbit pAb
Remarks/Alias	GLO2; GLX2; GLXII; HAGH1; HAGH
Species	Human
GeneID (human)	3029
GeneID	3029
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human HAGH (NP_001035517.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16775
Protein Weight	34kDa
Shipping	Ice bag